

Sleep Apnea Detector & Monitor



i-Sleep is a new device based on new prognostic index to screen and monitor sleep apnea conditions using the analysis of expired breath during sleep. i-Sleep is able to continuously assess respiratory functions with minimal patient effort in various sleep settings. i-Sleep comes with a screen displaying the novel sleep apnea index, end-tidal cardon dioxide (EtCO2) and respiration rate (RR). The patients only need to breathe comfortably during sleep to get the results. i-Sleep comes with a standard USB/SD card data storage powered by rechargeable batteries allowing data to be saved, monitored and analysed. The device is able to connect to smartphones via mobile apps where patients can independently monitor as well as manage their conditions.





Highlights

> Real-time Results

Measurements such as end-tidal carbon dioxide (EtCO2), respiratory rate (RR) and real-time capnogram were displayed within minutes.

> Portable and Lightweight

Designed to be easily managed and used, iBreath aims to deliver mobility and convenience to the user.

> Non-Invasive

i-Sleep offers a painless and non-invasive method to monitor your sleep apnea severity level.

> User-Friendly

Accompanied with simple steps that will effortlessly provide results within seconds. i-Sleep also caters to a wide range of ages from children to adults. Refer to user manual

> Cost-Effective and Reliable

Assembled with quality components that are cost-efficient, i-Sleep aims to deliver a quality performance at an affordable price. i-Sleep is capable of producing smooth EtCO2 waveform that can improve diagnosis and enhance its result reliability

> Data Storage

i-Sleep offers 16GB* worth of storage for data recording and saving.

> PC Based Via USB

i-Sleep's data can be accessed through a computer via a USB connection. With the data (in .csv format), users are free to replot, print and etc.

Performance

Physical Characteristics	Environmental	
Dimensionsw(140mm) xL(140m	nm)xH(89mm) Operating Temperature	0 - 40 degC
Weight Approx. 430g	Operating Humidity	10 - 90 % RH
Sample Rate	Storage Temperature	-40 - 70 degC
50 ml/minute	Storage Humidity	< 90 % RH
Measurement range	Pressure	At sea level
 0 - 114 mmHg. 0 - 15%. 0 - 15.2 kPa 	Battery TypeRechargeable	Lithium ion internal Battery
RR2 - 120 breath/minute (-	
0 - 40 mmHg: ±2 mmHg + 40 - 114 mmHg: ±5% of reading	1 Owel Bala	
	Battery Charger	AC 100V - 240V
	Output	5V 2A
RR ±1 BPM		

Specifications are subject to change without notice



- 1.Hjorth Activity
- 2. End-tidal Carbon Dioxide
- 3. Health Status
- 4. Exhaled CO2 Waveform
- 5. Respiratory Rate
- 6. Slope 3
- 7. Battery Percentage

Instructions

i-Sleep can produce results with this simple steps. (Applicable to device that has been completely set up)



2. Switch on the main switch and press the 'ON' button to start the device. Connect the nasal cannula to i-Sleep

 Properly wear the nasal cannula and position it comfortably with tape 3. View the waveform, HA, EtCO2, ICO2 and RR. Sleep normally at least 6 hours

4. Observe the health status that appear based on the analysed waveform.



For more elaborate instructions, please refer the user manual

Accessories



Nasal CO2 Sampling Cannula



Charging Adapter



Micro-SD Card with Adapter



Reference Pictures



i-Sleep w/ Nasal Cannula



Front



Right Side



Left Side

Product Availability

i-Sleep is currently available only in Malaysia.

Product Pricing

RMXXX

Purchase includes:
1x i-Sleep
1x Charging Adapter
1x Nasal CO2 Sampling Cannula
1x Micro SD Card (16GB) with Adapter

Contact Information

i-Respire Innovations Sdn. Bhd

No. 36, Jalan Pulai 22, Taman Pulai Utama, 81300, Skudai. Johor +607 - 5627290 info@i-respire.com

